


Radiometer in your hospital



Comprehensive acute care diagnostics

The Radiometer product portfolio

ABL800 FLEX analyzer



- Robust blood gas analyzer for high volume sites
- 18 STAT parameters from a single test
- Fully automated sample mixing and measurement

ABL90 FLEX analyzer



- Blood gas analyzer designed for ease of use at the POC
- Results in 35 sec from 65 µL samples
- Low maintenance and easy replacements

ABL80 FLEX analyzer series



- Blood gas analyzer ideal for low volume sites
- Small in size and fully operational on battery for high portability
- Cost-effective and flexible

Samplers



- Broad portfolio of syringes and capillaries
- Built-in features reduce risk of preanalytical errors and increase patient safety
- Lower risk of needle-stick injuries and no contact with patient blood for operator safety

Blood gas & related parameters

ANALYZER	TEST MENU																		
	pH	pCO ₂	pO ₂	ck ⁺	cNa ⁺	cCa ²⁺	cCl ⁻	cGlu	clac	cCrea	ctBili	sO ₂	ctHb	FO ₂ Hb	FCOHb	FMetHb	FHHb	FHbF	Hct
ABL90 FLEX	•	•	•	•	•	•	•	•	•		•	•	•	•	•	•	•	•	
ABL800 FLEX	•	•	•	•	•	•	•	•	•	•	•	•	•	•	•	•	•	•	•
ABL800 BASIC	•	•	•	•	•	•	•	•	•			•	•						
ABL80 FLEX CO-OX	•	•	•	•	•	•	•	•				•	•	•	•	•			
ABL80 FLEX	•	•	•	•	•	•	•	•											•
ABL80 FLEX BASIC version	•	•	•	•	•	•	•												•
ABL80 FLEX CO-OX OSM version												•	•	•	•	•	•		

AQT90 FLEX analyzer*



- Immunoassay analyzer
- High analytical performance
- Fast TAT for timely clinical decision making

IT systems



- One system for overview and control of POC devices
- Smart and scalable solution based on the newest Microsoft technology
- Intelligent compliance management and inspection readiness

TCM Series



- Continuous non-invasive transcutaneous monitoring of blood gas and patient status
- POC tcpO₂/tcpCO₂ monitoring of neonatal, pediatric and adult patients
- Diagnostic tool for patients suffering from vascular diseases

1st AUTOMATIC

- Simplified and automated blood sampling and data capture
- Reduced risk of preanalytical errors
- Secured sample integrity and operator safety
- Data accuracy: no patient-sample mix-ups, no lost billing opportunities



www.radiometer.com/1stautomatic

Immunoassay markers

TEST MENU \ ANALYZER	Troponin I	Troponin T	CKMB	Myoglobin	NT-proBNP	CRP	BhCG	D-dimer
AQT90 FLEX	•	•	•	•	•	•	•	•

TC parameters

TEST MENU \ MONITOR	tcpO ₂	tcpCO ₂	SpO ₂
TCM CombiM	•	•	
TCM TOSCA		•	•
TCM400	•		

Once implemented, we support your solution with training, compliance management and technical services. www.radiometer.com/services

* Not for sale in the US



ACUTE CARE TESTING