

# VITALOGRAPH PNEUMOTRAC

The Ideal PC-Based Spirometer for General Healthcare



**Vitalograph**<sup>®</sup>  
Your cardio-respiratory partner

Data you can rely on. People you can trust.

# VITALOGRAPH PNEUMOTRAC

The Ideal PC-Based Spirometer for General Healthcare and Occupational Health

The Vitalograph® Pneumotrac™ is the ideal PC-based spirometer for adult and paediatric use within general healthcare and occupational health. Pneumotrac is supplied complete with attractive hardware and software and connects directly to your laptop or PC.

Sophisticated and intuitive Vitalograph Spirotrac® Software provides the perfect environment for data collection, storage and trend analysis and includes the latest GLI predicted equations and Z-scores.

Integrated HL7 and GDT communication interfaces provide the platform for test data exchange with all common electronic medical record (EMR) systems




- Sophisticated PC-based spirometer
- Attractive, compact design
- Suitable for adults and paediatrics
- Durable measurement technology
- No need for expensive disposable sensors, turbines, spirottes or flow tubes
- No user calibration required
- Optimal hygiene via Vitalograph Bacterial Viral Filters (BVF™)
- Powerful Spirotrac Software with real-time curves, interpretation and data trending
- Choice of child friendly incentive animations to help improve patient cooperation
- Intuitive Spirotrac Software future-proofs your service through various PC expansion options including interfacing with your network and electronic medical record (EMR) systems
- Features latest GLI spirometry predicted equations and Z-scores

*“How fabulous your Pneumotrac is! They are the best results I have ever seen short of a rolling seal spirometer”*

Patrick F. Morgan, President | Morgan Scientific, Inc

Authorised Distributor



- ✓ Low Running Costs
- ✓ Optimal Hygiene

## Ordering Information

**77700 Vitalograph Pneumotrac Spirometer with Spirotrac Software**

## Accessories

- 28350 Vitalograph BVF Bacterial Viral Filters (50)
- 20303 Noseclips (200)
- 36020 Precision Syringe - 3 Litre

Vitalograph Ltd, UK. Maids Moreton, Buckingham, MK18 1SW, England. T: (01280) 827110 F: (01280) 823302 E: sales@vitalograph.co.uk www.vitalograph.co.uk  
 Vitalograph Ltd, Intl. Maids Moreton, Buckingham, MK18 1SW, England. T: +44 1280 827120 F: +44 1280 823302 E: sales@vitalograph.co.uk www.vitalograph.eu  
 Vitalograph Inc, 13310 West 99th Street Lenexa, Kansas, 66215, USA. Toll Free: 800 255 6626 T: (913) 730 3200 F: (913) 730 3232 E: vitcs@vitalograph.com www.vitalograph.com